


# Rheodyne® High Performance Fluidic Platforms

Raising the Bar in Precision Valve Technology



	TitanEX™	TitanHP™	TitanHT™
<b>Pressure Rating</b>	0-125 psi (0-9 bar) 	0-6000 psi (0-414 bar) 	0-15,000 psi (0-1034 bar) 
<b>Available Liquid Ends</b>	<ul style="list-style-type: none"> <li>▶ 2-position, 6-port</li> <li>▶ 6-position, 7-port</li> <li>▶ 10-position, 11-port</li> <li>▶ 2-position, 6-port, double 3-way</li> <li>▶ 2-position, 10-port</li> </ul>	<ul style="list-style-type: none"> <li>▶ 2-position, 6-port (analytical and nano)</li> <li>▶ 2-position, 10-port (analytical and nano)</li> <li>▶ vertical port, 2-position, 6-port</li> <li>▶ 6-position, 7-port</li> </ul>	<ul style="list-style-type: none"> <li>▶ 2-position, 6-port (analytical and semi-prep)</li> <li>▶ 2-position, 10-port</li> <li>▶ 6-position, 7-port selector</li> <li>▶ 3/7 injector with purge</li> <li>▶ 10-position, 11-port</li> <li>▶ 8-position, 4-port</li> </ul>
<b>Wetted Materials</b>	RPC-7*	DuraLife®, DuraLife II*, Stainless Steel w/PEEK™, PEEK	DuraLife, UltraLife™*, Stainless Steel w/PEEK, PEEK
<b>Port-to-Port Volume</b>	1.5-32.7 µL	0.28-0.50 µL	0.2-50.9 µL
<b>Flow Passage Diameters</b>	0.016-0.060 in 0.41-1.5 mm	0.004-0.012 in 0.10-0.30 mm	0.006-0.060 in 0.15-1.5 mm
<b>Operating Temperature</b>	0-60 °C	0-60 °C	0-60 °C
<b>Fittings</b>	Rheodyne TitanEX fittingless connections	10-32, M4	10-32, M4, 1/4-28
<b>Replacement Components</b>	Core assembly module (CAM) includes rotor, stator, and driver body	Rapid Replacement Pod™, RheBuild® Kit, replacement stators	Rapid Replacement Pod, RheBuild Kit, replacement stators, rotor seals

\* Rheodyne material combinations



	TitanEX™	TitanHP™	TitanHT™
<b>Suggested Applications</b>	Solvent selection, fraction collection, dispensation, process control/monitoring, IVD sample handling and preparation	High-pressure switching, injection, column selection, nano-scale applications, sample clean-up/enrichment, multidimensional chromatography	High-pressure switching, injection, column selection, sample clean-up/enrichment, multidimensional chromatography, and very-high-pressure (VHP) applications, prep to nano-scale applications
<b>Position Capability</b>	Multiposition with random access	Multiposition with random access	Multiposition with random access
<b>Optional Driver Board</b>	Yes	Yes	Yes
<b>Stop-to-Stop Actuation Speed</b>	280 ms	100–280 ms	100-280 ms
<b>Communication Protocols Available</b>	BCD, Level Logic Pulse, Dual Pulse, I <sup>2</sup> C, UART	BCD, Level Logic Pulse, Dual Pulse, I <sup>2</sup> C, UART	BCD, Level Logic Pulse, Dual Pulse, I <sup>2</sup> C, UART
<b>PCB Power Requirements</b>	24V @ 1 amp max	24V @ 1 amp max	Consult factory
<b>Motor Specifications</b>	Spark-free 7.5 degree stepper motor	Spark-free 7.5 degree stepper motor	Spark-free 1.8 degree stepper motor
<b>Overall Dimensions</b>	1.95 x 4.26 x 2.40 in. [49.5 x 108.2 x 61.0 mm]	1.95 x 4.87 x 2.40 in. [49.5 x 123.6 x 61.0 mm]	2.23 x 6.41 x 1.83 in. [56.6 x 162.8 x 46.6 mm]

DuraLife®, RheBuild®, Rapid Replacement Pod™, TitanEX™, TitanHP™, TitanHT™, and UltraLife™ are registered trademarks of IDEX Health & Science LLC. PEEK™ is a trademark of Victrex plc. © 2009 IDEX Health & Science LLC

R1-A4E:500-02/2009

